

ECX. CATION EXCHANGE POLYMER PHASE, STYRENE-DIVINYLBENZENE-BASED. IT ACTS AS A DUAL PHASE: REVERSED-PHASE AND CATION EXCHANGER.

PHASE	AMOUNT (mg)	VOLUME (ml)	PORUS (Å)	PARTICLE (µm)	PACK (U.)	ART. NO.
ECX	30	1	70	55	100	ECX030-01C
ECX	60	3	70	55	50	ECX060-03L
ECX	200	6	70	55	30	ECX200-06T
ECX	500	6	70	55	30	ECX500-06T

EAX. ANION EXCHANGE POLYMER PHASE, STYRENE-DIVINYLBENZENE-BASED. IT ACTS AS A DUAL PHASE: REVERSED-PHASE AND ANION EXCHANGER.

PHASE	AMOUNT (mg)	VOLUME (ml)	PORUS (Å)	PARTICLE (µm)	PACK (U.)	ART. NO.
EAX	30	1	70	55	100	EAX030-01C
EAX	200	6	70	55	30	EAX200-06T
EAX	500	6	70	55	30	EAX500-06T

SPE: BULK ADSORBENTS EXTRABOND®

BULK ADSORBENTS EXTRABOND®



ExtraBond® phases are available in bulk so you can prepare your own cartridges, or for extraction of solid samples by putting the packing directly in contact with the sample. Other phases are available and also larger quantities. Please contact helpdesk@scharlab.com.

PHASE	DESCRIPTION	PORUS (Å)	PACK (g)	ART. NO.
AL-N	Aluminium oxide neutral	60	100	000ALN100G
AL-N	Aluminium oxide neutral	60	1000	00ALN1000G
FL	Florisil	60	100	000FL0100G
FL	Florisil	60	1000	00FLO1000G
GCB	Graphitized carbon	-	100	CA03530100
GCB	Graphitized carbon	-	1000	CA03531000
NH2	Aminopropyl	60	100	000NH2100G
NH2	Aminopropyl	60	1000	00NH21000G
PSA	Ethyldiamine-N-propyl	60	100	KQ00230100
PSA	Ethyldiamine-N-propyl	60	1000	KQ00231000

PROBLEMS WITH THE METHOD OF EXTRACTION? IMPLEMENTING A NEW EXTRACTION PROCEDURE?



CONTACT SCHARLAB'S CHROMATOGRAPHY HELPDESK: helpdesk@scharlab.com

- WHAT IS THE MATRIX?
- WHAT ARE THE ANALYTES OF INTEREST?
- IS SPE BEING EMPLOYED FOR EXTRACTION?