

Aquagent® Titrant 2, free from pyridine, titrant-component for volumetric Karl Fischer titration**Identification**

Taric code: 3822 00 00

Applications

analytical chemistry, laboratory reagent.

Specifications

For use with: Aquagent® Solvent (AQ0002)

Contains methanol

Physical data

- Flash point: 11 °C
- Hygroscopic

Safety - GHS**Signal Word:** Danger**Hazard Statements:**

- H225: Highly flammable liquid and vapour.
- H331: Toxic if inhaled.
- H370: Causes damage to organs.

**Precautionary Statements:**

- P210: Keep away from heat / sparks / open flames / hot surfaces. - No smoking.
- P241: Use explosion-proof electrical / ventilating / lighting / equipment.
- P303+P361+P353: IF ON SKIN (or hair): Remove / Take off immediately all contaminated clothing. Rinse skin with water / shower.
- P370+P378: In case of fire: Use ... for extinction.
- P405: Store locked up.
- P501a: Dispose of contents / container in accordance with local / regional / national / international regulations.

Transport/storage

- ADR: 3 FT1 II • UN 1992 • FLAMMABLE LIQUID, TOXIC, N.O.S.(Aquagent® Titrant 2, titrant-component for Karl Fischer titration)
- IMDG: 3 II • UN 1992 • FLAMMABLE LIQUID, TOXIC, N.O.S.(Aquagent® Titrant 2, titrant-component for Karl Fischer titration)
- IATA/ICAO: 3 II • UN 1992 • FLAMMABLE LIQUID, TOXIC, N.O.S.(Aquagent® Titrant 2, titrant-component for Karl Fischer titration)
- PAX: 305
- CAO: 307
- 10°C - 30°C